

HIGH-SPEED CAMERA

Real-time Transmission Series



NEW!

1080P @ 1100FPS

HD Image Quality

72x72x60 (mm)

Compact Size

430g

Ultralight Camera

100ns

Ultra-short Exposure Time

Intelligent Image Triggered



Support Cross-frame PIV



Power Over CoaXPress



Automatic exposure
secondary exposure



Standardized CoaXPress-12
protocol. Real time transmission
speed > 4GB/S of image data



Support IRIG-B(DC)
time code input



SHS-201 is a new member of SinceVision intelligent high-speed camera family. The camera adopts advanced CXP-12 data interface, which can realize real-time writing speed > 2GB/s. With compact and sturdy body, outstanding performance and excellent heat dissipation, it can be handled calmly in various applications. With memory of computer, high-quality images and lossless data can be transmitted in real time, which is suitable for industrial scene monitoring, destructive experiments, system integration and other application scenarios.



Destructive Test



Precision manufacturing
monitoring



Microfluidics



Kinematics

Model: SHS-201

Max. Resolution: 1920 x 1080

Pixel size: 10μm

Max. Frame Rate of Full Frame: 1100FPS

Sensitivity: 30000 (M); 10000 (C)

Minimum Exposure Time: 100ns

Dynamic Range: 60dB

PIV Cross-frame Interval: ≥150ns

Color: Monochrome (M) / Color (C)

Trigger Mode: Manual / IO / Image

Customizable Storage: 2TB-30TB

Lens Mount: Optional E/C, adapt to F/EF

IO Signal: support TRIG/SYNC/IRIG-B

Data Interface: CXP-12, two channels

Protection Level: IP64

Bit Depth: 8bit/10bit

Power Supply Voltage: DC 12V-28V 12W

Working Temperature: 0-40°C, wide range version can be customized.


Camera Dimension: 60D x 72W x 72H

Camera Weight: 430g(E port), 470g(C port)

Frame Rate

Resolution	8bit	10bit
1920x1080	1120fps	850fps
1920x960	1260fps	1000fps
1920x720	1680fps	1300fps
1920x512	2350fps	1800fps
1920x256	4700fps	3700fps
1920x64	18900fps	15000fps

Minimal System Deployment


E-mail: info@sincevision.com

Website: www.sincevision.com

Quick Response Service Covering the Whole Country

Headquarters	: 5th Floor, Building 2, Chongwen Industrial Park, Nanshan Zhiyuan, Nanshan District, Shenzhen, China
East China Office	: Room 1305, Building 7, Xiangyu Liang'an Trade Center, No.1588, Chuangye Road, Kunshan, Jiangsu Province, China
North China Office	: Unit 922, Building 4, Times Fortune World, Courtyard 1, Hangfeng Road, Fengtai District, Beijing, China
Southwest China Office	: Room 604, Block B, Yingchuang International Building, No. 66, Chuangzhi South 1st Road, Pidu District, Chengdu, China
Northwest China Office	: Room 601, Chuangke Building, Cuihua Road, Yanta District, Xi 'an City, Shaanxi Province, China
Dongguan Office	: Room 406, Building F4, Tian'an Digital City, Nancheng District, Dongguan City, Guangdong Province, China
Hefei Office:	: Room 613, Tower C, Greenland Center, Hongcun Road, Baohe District, Hefei City, Anhui Province, China



SinceVision's
Linkedin Account



SinceVision's
YouTube Account